

AUTHORS:

Jin, Jianjian¹; Sabatino, Marianna¹; Somerville, Robert²; Wilson, John R.³; Dudley, Mark E.²; Stroncek, David F.¹; Rosenberg, Steven A.²

INSTITUTIONS:

1. Cell Processing Section, Department of Transfusion Medicine, Clinical Center, NIH, Bethesda, MD, United States.
2. Surgery Branch, National Cancer Institute, NIH, Bethesda, MD, United States.
3. Wilson Wolf Manufacturing Corporation, New Brighton, MN, United States.

CURRENT CATEGORY: Genetically Engineered Receptors and Adoptive Cell Therapies

PRESENTATION TYPE: Oral or Poster

ABSTRACT BODY: Adoptive cell therapy (ACT) of metastatic melanoma with autologous tumor infiltrating lymphocytes (TIL) is clinically effective, but TIL production can be challenging. Initial TIL culture and rapid expansion involves 4 to 6 weeks of culture in 24-well plates, tissue culture flasks and bags. Here we describe a simplified method for initial TIL culture and rapid expansion in gas-permeable flasks. TIL were initially cultured from tumor digests and fragments in 40 mL capacity flasks with a 10 cm² gas-permeable silicone bottom, G-Rex10. A TIL rapid expansion protocol (REP) was developed using 500 mL capacity flasks with a 100 cm² gas-permeable silicone bottom, G-Rex100. TIL growth was successfully initiated in G-Rex10 flasks from tumor digests from 13 of 14 patients and from tumor fragments in all 11 tumor samples tested. Enough TIL for REP could be obtained from one G-Rex10 flask (tumor digest) or 3 G-Rex10 flasks (tumor fragments) rather than two 24-well plates. TIL could then be expanded to 8-10 billion cells in a two-step REP which began by seeding 5 million TIL and 500 million irradiated peripheral blood mononuclear cells as feeder cells into a G-Rex100 flask. After 7 days of culture the cells were split and divided equally into 3 G-Rex100 flasks for the second step of TIL REP. This two-step G-Rex100 flask TIL REP was successful in 12 of 14 patients studied. To obtain the 30 to 60 billion cells used for patient treatment we seeded 6 G-Rex100 flasks with 5 million cells and later transferred the TIL into 18 G-Rex100 flasks. Large scale TIL REP in gas-permeable flasks requires approximately 9 to 10 liters of media, about 3 to 4 times less than other methods. In conclusion, TIL initiation and REP in gas-permeable G-Rex flasks require fewer total vessels, less media, less incubator space and less labor than initiation and REP in 24-well plates, tissue culture flasks and bags. TIL culture in G-Rex flasks will facilitate the production of TIL at the numbers required for patient treatment at most cell processing laboratories.